



The only clinically effective, non surgical, drug free treatment for glue ear

What is Otovent

Otovent autoinflation treatment is a clinically effective glue ear treatment designed to reduce the need for surgical intervention. It's convenient, increases chances of a shorter recovery time from glue ear and is the only effective non-invasive glue ear treatment for use during the 'watchful waiting' period.¹

Otovent may reduce ENT referrals for Glue Ear by up to 50%.


Otovent is listed on the drug tariff under Appliances: Autoinflation Device: Otovent
PIP Code 102889

NICE guidelines on Autoinflation

NICE recommends that GP's should offer autoinflation during the active observation period to children with OME who are likely to cooperate.²

Proven Success in Clinical Study³

Tympanometry⁴ after 2 weeks of Otovent treatment



About Glue Ear

Glue ear is a very common condition that affects 80% of children at some point during their childhood - 200,000 per year in the UK. Also known as secretory otitis media, otitis media with effusion or serous otitis media, more than 7 in 10 children have at least one episode of glue ear before they are 4 years old, and boys are more commonly affected than girls (source: Zeilhuis).

- 1 The time period between the first examination by the GP through to a grommet operation.
- 2 Quick reference guide. Issue date: February 2008. Surgical management of otitis media with effusion on children.
- 3 S-E Stangerup M.D., J. Sederberg-Olsen M.D., V. Balle M.D. Autoinflation as treatment of Secretory Otitis Media. Arch Otolaryngol Head Surg 1992; 118: 149-152.
- 4 The measurement of the outer and middle ear's ability to accept and conduct sound.